

RWD

Gas Filter Canister

Adsorb waste gases, guard health



R510-31-6



R510-31S-6



- ✓ Designed for the inhalation anesthesia scenarios of animals;
- ✓ Efficiently adsorbs anesthetic waste gases, including isoflurane, sevoflurane, enflurane, etc., preventing the direct discharge of anesthetics from polluting the environment;
- ✓ Single-use, with no risk of secondary pollution;
- ✓ Available in two sizes to meet the needs of various experimental settings;
- ✓ Adapted to various animal anesthesia environments such as mice, rats, rabbits, dogs and cats;
- ✓ Suitable for configuration methods such as passive type and active type;

Order Information:

Model	Product Description	Weight	Exhaust Gas Absorption Capacity	Remarks
R510-31-6	Gas Filter Canister, Large, pkg of 6pcs	800g	200g	Each can adsorb the waste gas produced by 200mL of isoflurane (two bottles)
R510-31S-6	Gas Filter Canister, Small, pkg of 6pcs	250g	50g	Each can adsorb the waste gas produced by 50mL of isoflurane (half a bottle)

- Paired with Gas Evacuation Apparatus, it ensures “0” gas leakage , making it safer.



Gas Evacuation Apparatus with Weight Monitoring (R546Pro)

- ✓ Provides an active exhaust gas recovery function and adopts an ultra-quiet motor;
- ✓ Supports multi-channel gas recovery, with continuously adjustable flow rate and real-time display on the LED screen;
- ✓ Features weighing function with audio-visual alarm for overloading: saturation warning at $\geq 990g$, overloading alarm at $\geq 1010g$;

One-Stop Professional Anesthesia Solution

► Multi-Function Anesthesia Solutions



Desktop Anesthesia Solutions



Mobile Anesthesia Solutions

Order Information:

Model	Product Description	Remarks
TAIJI-IE	Compact Small Animal Anesthesia Machine, Isoflurane, Easy Fill, 0-4L	Contains connection tubing for mask & induction chamber
TAIJI-SE	Compact Small Animal Anesthesia Machine, Sevoflurane, Easy Fill, 0-4L	
V105	Low Stress Anesthesia Induction Chamber, Mice and Rats (29cm*15cm*17cm)	Contains connection tubing for gas evacuation apparatus
V106	Low Stress Anesthesia Induction Chamber, Mice (16cm*10cm*11cm)	
R546Pro	Gas Evacuation Apparatus	
R510-29	Air Pump	Contains connection tubing for anesthesia machine
R510-T	Multi-Function Mobile Trolley	

► Isoflurane



Order Information:

Model	Product Description
R510-22-10	Isoflurane for Scientific Research Experiments, 100ml, pkg of 10pcs

One-Stop Professional Anesthesia Solution

► More Anesthesia Solutions




Small Animal Combination
Solutions of Ventilator to
Anesthesia System



Small Animal Stereotaxic
Anesthesia Solutions



Small Animal *in Vivo* Imaging
Anesthesia Solutions



RWD Life Science Co.,Ltd

Add: Building 7A (9th, 19th-20th Floors), Building 7D (9th Floor), Shenzhen International Innovation Valley, Dashi 1st Road, Nanshan District, Shenzhen, Guangdong, China
Enquiries: market@rwdstco.com

RWD Life Science Inc.

Add: 850 New Burton Road, Suite 201, Dover, DE 19904 (New Office)
Tel: (858)900-5879 E-mail: service@rwdls.com